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Hersteller / manufacturer: dent-e-con e.K.

Stand: 08/2018

Product: **TECHSIL 22 NF BLAU / BLUE**
Dubliersilikon / duplicating silicone

Instruction for use

Please read carefully and retain for future reference

Working
room temperature
at 23° C

Delivery form:
2x1 Kg bottles
2x5 Kg cans

2x6 Kg cans
2x25 Kg cans

Indications for use

Pourable, addition curing, two component silicone rubber that vulcanizes at room temperature.

Special features

- low viscosity
- easy to meter due to 1 : 1 ratio
- low hardness (22 Shore A)
- good mechanical properties
- very fast and non shrink cure at room temperature
- easy demolding

Important working hints

Silicones are high temperature sensitive. Please consider the recommended working and storage temperature.

Storage

TECHSIL is best stored between 5 °C and 30 °C in the tightly closed original container. The 'Best use before end' date of each batch is about 24 months after opening the container. Storage beyond the date specified on the label does not necessarily mean that the product is no longer usable. In this case however, the properties required for the intended use must be checked for quality assurance reasons.

Technical datas

Dosing:	1:1
Colour code:	red
Working time:	ca. 4 minutes
Setting time:	ca. 25 minutes
Hardness:	Shore A 22
Density:	1,08 g / cm ³
Viscosity:	ca. 8.000 mpA s
Linear dimension change:	< 0,2%
Tear resistance (DIN 53 504 S 1):	2,9 N/mm ²
Elongation at tear (DIN 53 504 S 1):	280 %
Continuing tear resistance (ASTM D 624 B):	3,8 N/mm

dent-e-con e.K.

Gartenstraße 19 | D-89173 Lonsee | T: +49 7336 9213369 | F: +49 7336 9213379

E: info@dent-e-con.de | W: www.dent-e-con.de